

Ear Canal Caps: General Information

1. Ear canal caps can be worn over the head, behind the head or under the chin.
2. Adjust the headband, where possible, to your approximate head size.
3. Place the headband in the desired position and readjust it, if necessary, for a snug comfortable fit. Each cap must block the ear canal opening completely.
4. Push the canal caps against the canal opening until sealed. A slight suction, when caps are removed, indicates a proper fit.
5. Do not bend the headband, since the tension is preset by the manufacturer. Wash with hand soap or a mild detergent.
6. **Suitable for wear only in steady-state noise levels below 95 dBA.**



DA Poster 40-501G, Jan 2005